# Elaiial®

Intra-Dermal Biotech HA

## DESCRIZIONE:

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## Intra-Dermal Biotech HA

### DESCRIPTION:


**ELAIAL** is a sterile, biodegradable and isotonic intradermal gel. **ELAIAL** consists of a medium molecular weight hyaluronic acid ( $1.5\text{-}2.0 \times 10^6$  Dalton), obtained from *Streptococcus equi* bacteria, formulated at a concentration of 20 mg/ml in a physiological buffer. Each box contains a **ELAIAL** syringe, two 30 G  $\frac{1}{2}$ " sterile, disposable needles reserved for the injection of **ELAIAL**, and a product leaflet. The package also contains two labels displaying the lot number. One of this labels should be affixed to the patient's medical file, while the other should be given to the patient to ensure traceability.

### COMPOSITION:

Sodium hyaluronate 20 mg/ml, sodium chloride, sodium phosphate monobase dihydrate, dibasic sodium phosphate dodecahydrate, WFI water. One syringe contains 2.0 ml of non-pyrogenic gel, sterilised using moist heat.

### INSTRUCTIONS FOR ASSEMBLING THE NEEDLE ON THE SYRINGE:

- Carefully unscrew the syringe cap, taking special care to avoid contact with the opening.
- Gently grasp the needle guard and fit the needle onto the Luer lock fitting, locking it tightly until a slight counterpressure is felt, to ensure an airtight seal and prevent the gel from leaking during administration.

### INDICATIONS:

**ELAIAL** is indicated in:

- Cases of scarring as a result of superficial skin lesions (e.g. acne and chickenpox scars).
- The physiological ageing process of the skin, whose effects include thickening of the stratum corneum (roughness and decreased skin radiance) and the alteration of the elastic fibres of the dermis.

### EXCLUSION CRITERIA:

**ELAIAL** must not be used in:

- Patients who tend to develop hypertrophic scarring.
- Patients with a history of autoimmune disease or who are undergoing immune therapy.
- Patients with known hypersensitivity to hyaluronic acid.
- Pregnant or breastfeeding women.
- Patients under 18 years of age.

Patients on anticoagulant or receiving platelet aggregation inhibitors should be treated with **ELAIAL** only after having consulted the attending physician.

Do not use **ELAIAL** in areas presenting inflammatory and/or infectious skin processes (e.g. acne, herpes, etc.). Do not use **ELAIAL** in association with laser therapy, chemical peeling or dermabrasion.

### PRECAUTIONS FOR USE:

**ELAIAL** is indicated exclusively for intradermal injections. Sensitive skin can be treated with creams or patches with local anaesthetic action. Please note that anaesthesia may cause local redness or hypersensitivity.

There is no clinical data (efficacy, tolerance) on the injection of **ELAIAL** into an area that has already been treated with another filling product. It is necessary to inform patients not to apply cosmetic products for 12 hours after the injection and to avoid prolonged exposure to sunlight and UV rays, nor to attend saunas or Turkish baths for a week after the injection. If the needle is obstructed, do not increase pressure on the plunger rod, but stop the injection and replace the needle. Do not implement injections into blood vessels. Do not inject **ELAIAL** into bones, tendons, ligaments or muscles. Do not inject **ELAIAL** into nevi (moles). After use, discard the syringe and any remaining product. Do not use product if the package is damaged.

### INTERACTIONS WITH OTHER AGENTS:

There are incompatibilities between sodium hyaluronate and quaternary ammonium compounds, such as benzalkonium chloride solutions. Therefore, avoid any contact between **ELAIAL** and these substances or medical-surgical instruments that have been in contact with these substances.

### SIDE EFFECTS:

Physicians must inform patients of potential side effects and/or incompatibilities associated with the implantation of this device, which may occur immediately after treatment or subsequently, and which include (the list is not exhaustive):

- Inflammatory reactions (e.g. redness, oedema, erythema, etc.) which may be associated with itching and pain on pressure after the injection. These reaction may last for a week.
- Hardening or nodules at the injection site, especially in cases of too superficial placement.
- Discolouration of the injection site.
- Insufficient effect if **ELAIAL** is injected incorrectly.

Patients should promptly inform their physician about any inflammatory reactions which persist for more than a week or any other developing side effects, so that they can be treated appropriately. Any other undesirable side effects associated with the injection of **ELAIAL** should be reported to the distributor and/or manufacturer.

### METHODS OF USE:

Inject **ELAIAL** into disinfected, healthy skin that is free of inflammation. The technique used is crucial to the success of the treatment. Therefore, the medical device should be used only in an adequately equipped clinic by physicians who have specific training in the injection technique. Before injection, thoroughly disinfect the area to be treated. Use the 30G  $\frac{1}{2}$ " needle supplied with the syringe and inject slowly by applying the appropriate injection technique. Assemble the syringe following the directions provided below. Inject **ELAIAL** at room temperature and under strict aseptic conditions. Inject **ELAIAL** intradermally with a linear or multi-picking (picotage) technique at a medium-deep level. The amount injected will depend on the wrinkle or skin depression to be treated and corrected. After injection, the physician may apply a gentle massage to evenly distribute the product. The use of **ELAIAL** is at the physician's discretion. **ELAIAL** should be used according to the needs of individual patients and the area to be treated.

### WARNINGS:

Before use, check the integrity of the syringe and the expiration date. Do not use needles other than those indicated. Do not reuse. As **ELAIAL** is a single-use product, quality and sterility are guaranteed only for an originally sealed syringe. Reuse of the product carries the potential risk of infection for patients or operators. Do not re-sterilize. After use, dispose of the needle in a suitable container.

**CAUTION:** the outside of the syringe is not sterile.

### STORAGE:

Store **ELAIAL** at 2-25°C (36-77°F) in a dry place in its original package. Protect from light, heat and extreme cold conditions. Keep out of the reach of children.

### FOR MEDICAL USE ONLY

### LAST REVISED ON:

12/2022

The Wave Innovation Group Srl  
Via M. Bellotti, 16 - 37139 Verona (Italy)  
www.thewaveinnovation.com - info@thewaveinnovation.com

### SYMBOLS

	Consult the package leaflet		Store between +2°C and +25°C
	Do not re-sterilize		Do not use if package is damaged
	Sterilized with moist heat		Do not reuse
	Keep in a dry place		Lot number
	Keep away from sunlight		Use by
	Sterile Fluid-Path Packaging		Product Reference Number
	Manufacturer		Date of production

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