# Produktil®

## Intra-Dermal Biotech HA

### INDICAZIONI:

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06/2022

The Wave Innovation Group S.r.l.  
Via M. Bellotti, 16 - 37139 Verona (Italy)  
www.thewaveinnovation.com - info@thewaveinnovation.com

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### INDICATIONS:

**PRODUKTIL** is indicated in the:

- Physiological ageing process of the skin, the effects of which include thickening of the stratum corneum (roughness and reduced skin radiance) and alteration of the elastic fibres of the dermis (wrinkles).
- Process of dermal tissue repair, in cases of scar sites due to superficial skin trauma (e.g. acne and chickenpox scars).

### DESCRIPTION:

**PRODUKTIL** is a sterile, biodegradable and isotonic intradermal gel. **PRODUKTIL** consists of a high molecular weight hyaluronic acid, obtained from *Streptococcus equi* bacteria and amino acids (L-proline, L-hydroxyproline, Glycine, L-lysine HCl), formulated at a concentration of 25 mg/ml in a physiological buffer. Each package contains two syringes of **PRODUKTIL** and a product leaflet. Inside the package, there are two sets of labels displaying the lot number. One of these labels should be affixed to the patient's medical file and the other should be given to the patient to ensure traceability.

### COMPOSITION:

Sodium hyaluronate 25 mg/ml, L-proline, L-hydroxyproline, Glycine, L-lysine HCL, Sodium phosphate monobase dihydrate, dibasic sodium phosphate dodecahydrate, WFI water.


### METHODS OF USE:

Inject **PRODUKTIL** into disinfected, healthy skin that is free of inflammation. The technique used is crucial to the success of the treatment. Therefore, the medical device should be used solely by physicians who have specific training in the injection technique. Before injection, thoroughly disinfect the area to be treated. The recommended needles to be used for this treatment are 27 to 30 G 1/2". After having disinfected the area, inject slowly by applying the appropriate injection technique. Assemble the syringe following the directions provided below. Inject **PRODUKTIL** at room temperature and under strict aseptic conditions. Inject **PRODUKTIL** intradermally with a linear or multi-picking (picoage) technique at a medium-deep level. The amount injected will depend on the wrinkle or skin depression to be treated and corrected. After injection, the physician may apply a gentle massage to evenly distribute the product. The use of **PRODUKTIL** is at the physician's discretion. **PRODUKTIL** should be used according to the needs of individual patients and the area to be treated.

### INSTRUCTIONS FOR ASSEMBLING THE NEEDLE ON THE SYRINGE:

A. Carefully unscrew the syringe cap, taking special care to avoid contact with the opening.

B. Gently grasp the needle guard and fit the needle onto the Luer-lock fitting, locking it tightly until a slight counterpressure is felt, to ensure an airtight seal and prevent the gel from leaking during administration.

### USE PRECAUTIONS:

**PRODUKTIL** is exclusively indicated for intradermal injections. Sensitive skin can be treated with creams or patches with local anaesthetic action. Please note that anaesthesia may cause local redness or hypersensitivity. There is no clinical data (efficacy, tolerance) on the injection of **PRODUKTIL** into an area that has already been treated with another filling product. It is necessary to inform patients not to apply cosmetic products for 12 hours after the injection and to avoid prolonged exposure to sunlight and UV rays, nor to attend saunas or Turkish baths for a week after the injection. If the needle is obstructed, do not increase pressure on the plunger rod, but stop the injection and replace the needle. Do not implement injections into blood vessels. Do not inject **PRODUKTIL** into bones, tendons, ligaments or muscles. Do not inject **PRODUKTIL** into nevi (moles). After use, discard the syringe and any remaining product. Do not use product if the package is damaged.

**CAUTION:** the outside of the syringe is not sterile.

### WARNINGS:

Before use, check the integrity of the syringe and the expiration date. Do not use needles other than those indicated. Do not reuse. As **PRODUKTIL** is a single-use product, quality and sterility are guaranteed only for an originally sealed syringe. Reuse of the product carries the potential risk of infection for patients or operators. Do not re-sterilize. After use, dispose of the needle in a suitable container.

### EXCLUSION CRITERIA:

**PRODUKTIL** must not be used in:

- Patients who tend to develop hypertrophic scarring.
- Patients with a history of autoimmune disease or who are undergoing immune therapy.
- Patients with known hypersensitivity to hyaluronic acid.
- Pregnant or breastfeeding women.
- Patients under 18 years of age.

Patients on anticoagulant or receiving platelet aggregation inhibitors should be treated with **PRODUKTIL** only after having consulted the attending physician. Do not use **PRODUKTIL** in areas presenting inflammatory and/or infectious skin processes (ex. acne, herpes, etc.). Do not use **PRODUKTIL** in association with laser therapy, chemical peeling or dermabrasion.

### INTERACTION WITH OTHER AGENTS:

There are incompatibilities between sodium hyaluronate and quaternary ammonium compounds, such as benzalkonium chloride solutions. Therefore, avoid any contact between **PRODUKTIL** and these substances or medical-surgical instruments that have been in contact with these substances.

### SIDE EFFECTS:

Physicians must inform patients of potential side effects and/or incompatibilities associated with the implantation of this device, which may occur immediately after treatment or subsequently, and which include (the list is not exhaustive):

- Inflammatory reactions (e.g. redness, edema, erythema, etc.) which may be associated with itching and pain on pressure after the injection. These reaction may last for a week.
- Hardening or nodules at the injection site, especially in cases of too superficial placement.
- Discolouration of the injection site.
- Insufficient effect if **PRODUKTIL** is injected incorrectly.

Patients should promptly inform their physician about any inflammatory reactions which persist for more than a week or any other developing side effects, so that they can be treated appropriately. Any other undesirable side effects associated with the injection of **PRODUKTIL** should be reported to the distributor and/or manufacturer.

### STORAGE:

Store **PRODUKTIL** at 2-25°C (36-77°F) in a dry place in its original package. Protect from light, heat and extreme cold conditions. Keep out of the reach of children.

### PACKAGE CONTENT:

2 pre-filled syringes, each containing 2ml of non-pyrogenic gel, sterilised using moist heat.

### FOR MEDICAL USE ONLY

### LAST REVISION ON:

06/2022

### SYMBOLS

	Consult the package leaflet		Store between +2°C and +25°C
	Do not re-sterilize		Do not use if package is damaged
	Sterilized with moist heat		Do not reuse
	Keep in a dry place		Lot number
	Keep away from sunlight		Use by
	Sterile Fluid-Path Packaging		Product Reference Number
	Manufacturer		Date of production

The Wave Innovation Group Srl  
Via M. Bellotti, 16 - 37139 Verona (Italy)  
www.thewaveinnovation.com - info@thewaveinnovation.com

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