Produktil®

Intra-Dermal Biotech H A

INDICAZIONI:

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The Wave Innovation Group srls
Via Barucchi, 47/A - 37139 Verona (Italy)

www.thewaveinnovation.com - info@thewaveinnovation.com

0373

Produktil®

Intra-Dermal Biotech H A

INDICATIONS:

PRODUKTIL is indicated:

- in the physiological process of ageing of the skin, the effects of which include inspissation of the horny layer (roughness and reduced radiance of the skin) and changes in the elastic fibres of the dermis (wrinkles).
- in the dermal tissue repair process, in cases of scar results following superficial cutaneous trauma (e.g. acne and chickenpox scars).

DESCRIPTION:

PRODUKTIL is a sterile, biodegradable, isotonic, intradermal gel. **PRODUKTIL** consists of medium chain hyaluronic acid, obtained from *Streptococcus equi* bacteria, and amino acids (L-proline, L-hydroxyproline, Glycine, L-lysine HCl), formulated to a concentration of 25 mg/ml in a physiological buffer. Each box contains two syringes of **PRODUKTIL** and a product leaflet. A set of two couples of labels showing the batch number is contained in the box. One of these labels should be attached to the patient's file and the other should be given to the patient to ensure traceability.

COMPOSITION:

Sodium hyaluronate 25 mg/ml, L-proline, L-hydroxyproline, Glycine, L-lysine HCL, sodium dihydrogen phosphate dihydrate, disodium hydrogen phosphate dodecahydrate, WFI grade water.


METHODS OF USE:

PRODUKTIL must be injected into non-inflamed, disinfected, healthy skin. The technique used is essential for the success of the treatment. Therefore this device must only be used by doctors who have received specific training on the injection technique. The area to be treated should be disinfected thoroughly prior to the injection. We recommend the use of needles from 27 to 30G ½", so inject slowly by applying the appropriate injection technique. Inject **PRODUKTIL** at room temperature and with strict asepsis conditions. Inject **PRODUKTIL** intradermally with a linear technique or with picotage at a medium deep level. The amount injected will depend on the wrinkle or depression which is to be corrected. After the injection, doctors may apply a light massage in order to distribute the product uniformly. The use of **PRODUKTIL** is at the discretion of the doctor. **PRODUKTIL** must be used according to the needs of individual patients and the area to be treated.

INSTRUCTIONS FOR ASSEMBLY OF THE SYRINGE NEEDLE:

A. Carefully unscrew the cap of the tip of the syringe, being particularly careful to avoid contact with the opening.

B. Gently grip the needle guard and mount the needle on the luer-lock mount, screwing it tight until a slight counter-pressure is felt in order to ensure an airtight grip and prevent leakage of the liquid during administration.

PRECAUTIONS FOR USE:

PRODUKTIL is only indicated for intradermal injections. Sensitive skin may be pretreated using a local anaesthetic patch or cream. Please note that any anaesthesia may cause redness or local hypersensitivity. There are no available clinical data (efficiency, tolerance) about injecting **PRODUKTIL** into an area which has already been treated with another filling product. Patients should be advised not to apply any make-up for 12 hours after the injection and to avoid prolonged exposure to sunlight and UV or using saunas or Turkish baths for one week after the injection. If the needle is blocked, do not increase the pressure on the plunger rod but stop the injection and replace the needle. Do not inject into blood vessels. Do not use **PRODUKTIL** in bones, tendons, ligaments or muscles. Do not inject **PRODUKTIL** into naevi. Discard the syringe and remaining product after use. Do not use if packaging has been damaged.

WARNINGS:

Verify the integrity of the syringe and the expiry date before use. Do not use needles other than those listed. Do not re-use; quality and sterility can only be guaranteed for an originally closed syringe. The re-use of the product creates a potential infection risk for patients or users. Do not re-sterilise. After use the needle must be disposed of in a suitable container.

EXCLUSION CRITERIA:

PRODUKTIL must not be used in:

- patients who tend to develop hypertrophic scarring,
- patients with a history of autoimmune disease or who are receiving immune therapy,
- patients who are known to be hypersensitive to hyaluronic acid,
- pregnant or breastfeeding women,
- patients under 18 years of age.

Anticoagulated patients or patients receiving platelet aggregation inhibitors should not be treated with **PRODUKTIL** without consulting their doctors. **PRODUKTIL** must not be used in areas presenting cutaneous, inflammatory and/or infectious processes (e.g. acne, herpes ...). **PRODUKTIL** must not be used in association with laser therapy, chemical peeling or dermabrasion.

INTERACTIONS WITH OTHER AGENTS:

There are incompatibilities between sodium hyaluronate and quaternary ammonium compounds such as benzalkonium chloride solutions. Therefore **PRODUKTIL** should never be placed in contact with these substances or with medical-surgical instruments that have been in contact with these substances.

UNDESIRABLE SIDE EFFECTS:

Physicians must inform the patient that there are potential side effects and/or incompatibilities associated with the implantation of this device, which may occur immediately or may be delayed. These include (non-exhaustive list):

- Inflammatory reactions (e.g. redness, oedema, erythema ...) which may be associated with itching and pain on pressure after the injection. These reactions may last for a week.
- Inclusions or nodules at the injection site especially in cases of too superficial placement.
- Discolouration of the injection site.
- Poor effect if **PRODUKTIL** is injected incorrectly.

Patients must inform their physician as soon as possible about any inflammatory reactions which persist for more than one week or any other secondary effect which develops. The physician should treat these as appropriate. Any other undesirable side effects associated with the injection of **PRODUKTIL** must be reported to the distributor and/or to the manufacturer.

STORAGE:

Store **PRODUKTIL** at 2–25°C (36–77°F) in a dry place in the original box. Protect from light, heat and frost and handle with care.

CONTENTS OF THE PACK:

2 pre-filled syringes each containing 2 ml of non pyrogenic gel, sterilised using moist heat.

ONLY FOR MEDICAL USE

LAST REVISED:

01/2020

SYMBOLS

 Consult instructions for use

 Do not re-sterilise

 Sterilised using steam

 Keep dry

 Keep away from sunlight

 Sterile fluid path

 Store between +2°C and +25°C

 Do not use if package is damaged

 Do not reuse

 Manufacturer

 Batch number

 Use by

The Wave Innovation Group srls
Via Barucchi, 47/A - 37139 Verona (Italy)

www.thewaveinnovation.com - info@thewaveinnovation.com

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