


Profiler


Hyaluronic Dermal Filler

THE WAVE INNOVATION GROUP presents an innovative fillers based on bio fermentative hyaluronic acid sourced no animal origin. Profiller is made thanks the latest studies and biotechnological systems that ensuere maximum security and reliability combing excellent performance and quality results



PROFILLER, is made with a reticulate highly linked to a pharmacological hyaluronic acid to give more stability and durability.


It is a hypoallergenic product formulated at a concentration of 20 mg / ml, crosslinked with the lowest percentage of BDDE and endotoxins (<0.5 EU / G) compared to the best filling in the market. In fact, you do not need allergy tests before treatment.



PROFILLER: It is a soft, homogeneous, biodegradable, sterile, viscoelastic, clear, transparent and isotonic gel injection.

PROFILLER is a single phase filling with the best relation between time, stability and comfort for the patient during the application.


Developed to fill facial wrinkles, furrows, small depressed areas, with great utility to increase the volume of the cheekbones and lips, with reabsorption times ranging from 8 months of application onwards.



Switch cross-linking technology:

Profiller employs a sophisticated and advanced purification process, alignment and crosslinking Switch of hyaluronic acid, this procedure allows to obtain a homogenous thick gel because there is no passage of other particles through the alignment of the molecular chains.

Our studies shows a duration of 6-8 months into wrinkles for Profiller Easy and 7-10 months for Profiller Lips and Visage.



| PROFILLER | EASY | VISAGE | LIPS |
|------------------------|--|---|-------------------|
| Indicated | Peri ocular área, Peri bocal area (barcode) Micro wrinkles | Tunneling wrinkles, Filling of the chin, Nasogenian wrinkles, Cheekbones | Lips Cupid bow |
| Content | 1 ml | 1 ml | 1 ml |
| Concentration HA | 20 mg/ml | 20 mg/ml | 20 mg/ml |
| Reticulate agent | BDDE | BDDE | BDDE |
| Degree of Crosslinking | ***** | ***** | ***** |
| Needles | 2/30 G | 2/27 G | 2/27 G |
| Duration | 6-8months | 7-10months | 7-10 months |

ONE DERMAL FILLER FOR EACH SOLUTION



Smooth brow furrows, crow's feet, frown lines
Profiller Visage



Smooth under the eyes/ tear troughs
Profiller Easy



Lift & enhance cheeks
Profiller Visage



Soften nasolabial folds
Profiller Visage

Lips & smooth vertical lips lines
Profiller Lips



Smooth marionettes lines
Profiller Visage




Smooth chin wrinkles
Profiller Visage



PROTOCOL / DOSAGE

Making a session in the
area of interest



Profiller

Hyaluronic Dermal Filler



www.thewaveinnovation.com